





# Up to now ( $\leq$ VB17A) – File Structure

## meas.dat




BufferHeader  
(variable length)





## Up to now ( $\leq$ VB17A) – One Scan

### ■ Measurement header + data




\\n4\pkg\MrServers\MrMeasSrv\SeqIF\MDH\mdh.h



# Motivation for Multi-RAID


- Remove redundant MDH info (receiver performance)
  - Extend MDH with additional application parameters
- 
- Before VD11A only a single measurement data can be copied by TWIX
    - Problem: If PreScan-Normalize or NoiseDecorrelation is used also these pre-measurement data have to be copied to run ICE simulation.
    - Solution: New Multi-RAID file in VD11A (selectable in TWIX) collects automatically all dependent measUIDs needed for ICE simulation.





## Multi-RAID – file structure

### ■ meas.dat (VD11A)




- The last measurement contains the „real“ image scan data, first one's are adjustments, pre-scans, noise-scans
- The last measurement depends on the preceding measurements
- Use offset to determine the position of a measurement in file
- There are always 64 MrParcRaidFileEntry structs, independent of the actual number of measurements
- After the MultiRaidFileHeader and after every measurement there is a 512 byte alignment filled with „0“



## Multi-RAID – file sub structure

### ■ VD11A:




unchanged structure    new structure



## Multi-RAID – One Scan

- VD11A: MDH = ScanHeader + ChannelHeader




VB17A: FreePara → use IceProgramPara in future  
 FreePara[0]..[3] ≡ IceProgramPara[4]..[7]  
 VD11A: ReservedPara, must not be used by any application

\\n4\pkg\MrServers\MrMeasSrv\SeqIF\MDH\mdh.h



## Multi-RAID – Pixel data


- VB17A  $\equiv$  VD11A:



4 Byte | 4 Byte  
Re{z} Im{z}



## RAID file – Summary: Old vs. New







## Up to now ( $\leq$ VB17A) – MDH

\n4\pkg\MrServers\MrMeasSrv\SeqIF\MDH\mdh.h

<b>DMALength</b>	DMA length [bytes] (bit 0..24 of 4 Byte)
<b>MeasUID</b>	measurement user ID (4 Byte)
<b>ScanCounter</b>	scan counter [1...] (4 Byte)
<b>TimeStamp</b>	time stamp [2.5 ms ticks since 00:00] (4 Byte)
<b>PMUTimeStamp</b>	PMU time stamp [2.5 ms ticks since last trigger] (4 Byte)
<b>EvalInfoMask[1..2]</b>	evaluation info mask field (2 · 4 Byte)
<b>SamplesInScan</b>	# of samples acquired in scan(2 Byte)
<b>UsedChannel</b>	# of channels used in scan(2 Byte)
<b>sLoopCounter</b>	loop counters [1..7] (28 Byte)
<b>sCutOffData</b>	cut-off values(4 Byte)
<b>KSpaceCentreColumn</b>	centre of echo(2 Byte)
<b>CoilSelect</b>	Coil select (2 Byte)
<b>ReadOutOffcentre</b>	ReadOut offcenter value (4 Byte)
<b>TimeSinceLastRF</b>	Sequence time stamp since last RF pulse (4 Byte)
<b>KSpaceCentreLineNo</b>	number of K-space centre line(2 Byte)
<b>KSpaceCentrePartitionNo</b>	number of K-space centre partition(2 Byte)
<b>IceProgramPara</b>	free parameter for IceProgram (8 Byte)
<b>FreePara</b>	free parameter (8 Byte)
<b>SliceData</b>	Slice Data (28 Byte)
<b>ChannelId</b>	channel Id (2 Byte)
<b>PTABPosNeg</b>	negative, absolute PTAB position in [0.1 mm] (2 Byte)

$\Sigma$  = 128 Byte



## VD11A – ScanHeader

\\n4\pkg\MrServers\MrMeasSrv\SeqIF\MDH\mdh.h

<b>DMALength</b>	DMA length [bytes] (bit 0..24 of 4 Byte)
<b>MeasUID</b>	measurement user ID (4 Byte)
<b>ScanCounter</b>	scan counter [1...] (4 Byte)
<b>TimeStamp</b>	time stamp [2.5 ms ticks since 00:00] (4 Byte)
<b>PMUTimeStamp</b>	PMU time stamp [2.5 ms ticks since last trigger] (4 Byte)
<b>SystemType</b>	System type (2 Byte)
<b>PTABPosDelay</b>	PTAB delay (2 Byte)
<b>PTABPosX</b>	absolute PTAB position in [0.1 mm] (4 Byte)
<b>PTABPosY</b>	absolute PTAB position in [0.1 mm] (4 Byte)
<b>PTABPosZ</b>	absolute PTAB position in [0.1 mm] (4 Byte)
<b>Reserved1</b>	reserved for future hardware signals (4 Byte)
<b>EvalInfoMask[1..2]</b>	evaluation info mask field (2 · 4 Byte)
...<like in VB17A>	...(50 Byte)
<b>KSpaceCentrePartitionNo</b>	number of K-space centre partition(2 Byte)
<b>SliceData</b>	Slice Data (28 Byte)
<b>IceProgramPara</b> ▶	free parameter for IceProgram (48 Byte)
<b>ReservedPara</b>	unused internal reserved parameter (8 Byte)
<b>ApplicationCounter</b>	(2 Byte)
<b>ApplicationMask</b>	(2 Byte)
<b>CRC</b>	CRC 32 checksum (4 Byte)

Σ = 192 Byte



## VD11A – ChannelHeader

\\n4\pkg\MrServers\MrMeasSrv\SeqIF\MDH\mdh.h

<b>TypeAndChannelLength</b>	DMA length [bytes] (bit 8..31 of 4 Byte)
<b>MeasUID</b>	measurement user ID (4 Byte)
<b>ScanCounter</b>	scan counter [1...] (4 Byte)
<b>Reserved1</b>	reserved ( 4 Byte)
<b>SequenceTime</b>	sequence readout starting time (4 Byte)
<b>Unused2</b>	unused (4 Byte)
<b>ChannelId</b>	channel Id (2 Byte)
<b>Unused3</b>	unused (2 Byte)
<b>CRC</b>	CRC 32 checksum of channel header (4 Byte)

Σ = 32 Byte



## VD11A – MultiRaidFileHeader

\n4\pkg\MrServers\MrVista\include\Parc\MrParcRaid.h

### MrParcRaidFileHeader

#### HeaderSize

Used for MDH file identification (4 Byte)

< 32: VD11A MDH file

≥ 32: Old pre VD11A MDH file

Values:

	MR_PARC_RAID_ALLDATA	0	normal VD11A file (*.dat)	} VD line
	MR_PARC_RAID_MDHONLY	1	compact file (meas'data removed) (*.mdh)	
	MR_PARC_RAID_HDONLY	2	file only with Multi-RAID and buffer header	
	MR_PARC_RAID_LEGACY_THR	32	pre VD11A file without buffer header (no RAID)	VA line
	> 32		pre VD11A file with buffer header	VB line

#### MeasCount

Number of measurements in file (4 Byte)

Σ = 8 Byte

### MrParcRaidFileEntry

#### Measurement UID

Measurement UID (4 Byte)

#### File ID

File ID from RAID (4 Byte)

#### Offset

Offset to measurement (and header) from start of file (8 Byte)

#### Length

Length of measurement and header data (8 Byte)

#### PatName

Patient name (64 Byte)

#### ProtName

Protocol name (64 Byte)

Raid mode - Twix	
FID	MeasID
MID107_gre_FID107	107 107
MID106_gre_FID106	106 106
MID105_gre_FID105	105 105

Σ = 152 Byte